

## GE Fight Bac™ Disinfectant Wipes (Canada)



### Disinfectant Wipes

Disinfectant with the environment in mind. GE Fight Bac™ Wipes kills over 99.9% of bacteria and viruses in 5-minutes of contact time. Formulated with a botanical active ingredient, citric acid, and free from harsh chemicals, this product provides superior disinfection

[Order GE Fight Bac Wipes Today!](#)

[Search Betco on Health Canada's list of disinfectants with evidence for use against COVID-19.](#)

### Features

Meets Health Canada's criteria for use against SARS-CoV-2, the virus that causes COVID-19  
Hospital grade efficacy  
Citric Acid active ingredient is non-irritating and no harsh fumes for the user  
Environmentally friendly

### Available In

4 - 500 Wipes Bucket 5.5" x 7" sheet Item #392F1-07

### Directions for Use

**In the presence of 5% blood serum and a contact time of 5 minutes, is effective against:**

**BACTERICIDAL:** Pseudomonas aeruginosa Staphylococcus aureus Salmonella enterica, Listeria monocytogenes, Escherichia coli O157:H7, Staphylococcus aureus - MRSA, Staphylococcus epidermis - MRSE, Enterococcus faecalis Vancomycin Resistant (VRE)

**VIRUCIDAL:** Respiratory Syncytial Virus, Rotavirus, Herpes Simplex Type 1 Virus, Herpes Simplex Type 2 Virus, Vaccinia Virus, Influenza A Virus, Adenovirus Type 2, Feline Calicivirus (surrogate for Norovirus), Human Hepatitis B Virus (HBV)

**FUNGICIDAL:** Trichophyton interdigitale.

&#160;

**In the presence of 5% blood serum and a contact time of 10 minutes, is effective against:**

**BACTERICIDAL:** Acinetobacter baumannii

**VIRUCIDAL:** Canine Parvovirus, Poliovirus Type 1

This product is a broad-spectrum virucidal hard surface disinfectant that is expected to inactivate the SARS-CoV-2 (the virus that causes COVID-19)

**In the presence of 5% blood serum and a contact time of 1 minute, is effective against:** Human Immunodeficiency Virus Type 1

**MYCOBACTERICIDAL:** Effective against Mycobacterium bovis in 5 minutes at 20&#176;C / 68&#176;F on previously cleaned hard non-porous inanimate surfaces.

**SANITIZER:** Staphylococcus aureus and Klebsiella pneumonia in 1 minute (60 seconds).

**AREAS OF USE:** Hospitals, Kitchens, Bathrooms, Health Care Facilities, Medical Offices, Day Care Centers, Schools and Nursing Homes

**SURFACES:** Floors, Fixtures, Walls, Countertops, Bathrooms, Computers, Metal, Stainless Steel, Plastic, Vinyl and Glass

### PRECAUTIONARY STATEMENTS

#### HAZARD TO HUMANS AND DOMESTIC ANIMALS

**CAUTION:** Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet.

### FIRST AID

**IF IN EYES:** Hold eyelids open and flush with a steady, gentle stream of water for 15 minutes. Remove contact lenses, if present, after the first five minutes, then continue rinsing. Call a Poison Control Center or doctor for further treatment advice.

**DIRECTIONS FOR USE:** Product is ready to use. Preclean areas with gross filth and heavy soil prior to disinfection.

**DISPENSER DIRECTIONS:** To start feed, remove large cover. Pull up the towelette corner from the center of the towelette roll, twist it into a point and thread it through the hole in

the container cover. Pull through about one inch. Replace cover. Pull out first towelette and snap off at a 90° angle. Remaining towelettes feed automatically ready for their next use. When through using, keep small center cap closed to prevent moisture loss.

**TO CLEAN/DEODORIZE:** Unfold premoistened cloth and wipe area to be cleaned for 30 seconds or until clean. Spot check to test for compatibility

**TO SANITIZE:** Unfold premoistened cloth and wipe hard non-porous area or item to be sanitized. Wet all surfaces thoroughly. Let stand for 60 seconds (one minute). Wipe or let air dry. No rinsing or wiping is required, except on direct food contact surfaces which require a potable water rinse after treatment.

**TO DISINFECT:** Unfold premoistened cloth and wipe hard non-porous area or item to be disinfected. Wet all surfaces thoroughly. The surface must remain wet with GE Fight Bac Wipes for at least 5 or 10 minutes. Wipe or let air dry. No rinsing or wiping is required, except on direct food contact surfaces which require a potable water rinse after treatment. Toys: Use only on hard non-porous surfaces and rinse after treatment. **FUNGICIDAL DIRECTIONS:** This product is effective against Trichophyton interdigitale (athlete's foot fungus and a cause of ringworm) on inanimate surfaces in locker rooms, dressing rooms, exercise facilities, and shower and bath areas. Follow disinfection directions.

**Rinse toys and food contact surfaces with potable water.**

**DECONTAMINATION OF CRITICAL OR SEMI-CRITICAL DEVICES PRIOR TO TERMINAL STERILIZATION/HIGH LEVEL DISINFECTION:** This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or, (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

**GE FIGHT BAC WIPES KILLS HIV and HBV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in health care settings (Hospitals, Nursing Homes) or other settings in which there is an expected likelihood of soiling of hard, non-porous surfaces/objects with blood or body fluids and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of human immunodeficiency virus Type 1 (HIV-1) (associated with AIDS) or Human Hepatitis B Virus (HBV).

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATING SURFACES/OBJECTS SOILED WITH BLOOD AND BODY FLUIDS WHICH MAY CONTAIN HIV-1 (AIDS virus) or Human Hepatitis B Virus (HBV).**

**PERSONAL PROTECTION:** Personnel that clean items soiled with blood or body fluids must wear appropriate personal protective equipment (i.e., barrier protection), such as disposable gloves, gowns, masks or eye coverings.

**CLEANING PROCEDURE:** Target surfaces must be cleaned prior to disinfection, as the presence of heavy soil on a surface may reduce the intended disinfectant efficacy of the product.

**CONTACT TIME:** Leave precleaned surfaces wet for 1 minute to mitigate HIV-1 virus or 5 minutes to mitigate Human Hepatitis B Virus (HBV).

**DISPOSAL OF INFECTIOUS MATERIAL:** Any materials used in the cleaning process that may contain blood or body fluids are to be disposed of immediately in accordance with local

regulations for infectious materials disposal.

**STORAGE:** Do not contaminate water, food or feed by storage and disposal. Do not leave containers open during storage. No special handling is required.

**DISPOSAL: Towelette:** Do not reuse towelette. Dispose of used towelette in trash. Dispenser: Nonrefillable container. Do not reuse or refill this container. Offer for recycling if available.

**Tech Specs**

Fragrance	Fresh
pH	2-2.3
Wipe color	White
DIN	02513943
Warranty	1 year