

Symplicity™ PRO G Healthcare Prewash



Healthcare Prewash

Ready-to-Use prewash spotter specifically formulated to remove many of the stubborn stains commonly found in the longterm and acute healthcare market including blood, Granulex, Betadine and other medicinal based stains. When used as directed and in accordance with Betco's Best Laundry Practices wall chart, PRO G will save time energy labor and reduce chemical consumption and produce one wash clean results. Part of Betco's LaundryPro Stain Treatment Center

Features

Convenient packaging – Ready to Use
Reduces Rewash - Labor Saving
Designed for stains found the Healthcare environments
Fast Acting – Labor Savings

Available In

6 - 32 oz. Oval Bottles Item #48816-00

Directions for Use

Spray Pro G directly onto stain. Agitate with brush or rub vigorously. If possible, wait 5 to 10 minutes before beginning wash cycle. Always follow label directions.

IMPORTANT: Before using, test for dye stability by cleaning a small inconspicuous section of fabric.

Tech Specs

| | |
|-----------|-------------|
| Color | Colorless |
| Fragrance | Country Air |
| Foaming | Slight |
| Warranty | 2 years |