

## Crete Rx™ CRP



### Crete Rx CRP

This dual action counter-rotating plate attaches directly to standard floor machines, enabling the user to increase the performance and productivity when working on concrete floor surfaces with the Crete Rx Honing and Polishing prescription pathway.

### Features

Added 26 Lbs of head pressure significantly reduces time and labor  
Allows Floor Machines to better remove stains and etches in Concrete Floors & produce higher depth of gloss  
Durable Aluminum & Iron Construction is long lasting, for years of Crete Rx projects

### Available In

**Crete Rx™ CRP Item #E84143-00**

### Directions for Use

Attach Crete Rx CRP to any standard 17" or 20" floor machine. For stain and etch remediation, make 4 to 5 passes over the stained/etched area at a speed of 3 to 4 feet/minute with the Crete Rx Honing Tool. During the Polish step, make 4 to 5 passes over the stained/etched area at a speed of 3 to 4 feet/minute with the Polishing Tool

### Tech Specs

Disc Width	16"
Satellite Width	3"
Weight	25 Lbs
Construction	Iron & Aluminum